

BMP®21-LAB Label Printer



bradyID.com.au/bmp21lab

Where
Pens Fail.
Printers
Succeed.

When it comes to in-the-lab labelling, illegible or smeared writing is not an option. With the powerful and affordable BMP®21-LAB printer on hand, you can rely on crisp, long-lasting identification. Once you experience the built-in symbols, durable labels and vial/tube auto-sizing capabilities of the BMP21-LAB printer, you'll never use a pen again.

Labels Tested for Labs

- Labels and inks built to withstand harsh lab environments, including autoclaves, freezers, liquid nitrogen, and isopropyl alcohol and methyl ethyl ketone chemicals
- Unlike handwritten identification, printed labels reduce errors with easy-to-read text that fits more information in a small space
- Labels engineered for the lab means you never lose your data due to a failed label
- Technical Data Sheets (TDS) available with the specific properties and test results of each material, available at bradyID.com.au/techdata

Convenient Printer Features

- Vial label shortcut key automatically cuts labels to fit to vial and tube sizes
- Greek symbols and time and date stamps for enhanced sample tracking
- Built-in serialization and bar code symbologies provide more information in a smaller footprint
- Drop-lock-and-print cartridges for quick printing with no setup required
- Easy to use and easy to train others to use
- Variety of power options (use AA batteries, rechargeable battery or AC power supply — all sold separately)

Easy to Own and Maintain

- Label supplies and printer available through top lab distributors
- Low price point allows each technician or lab station to have a printer
- Technical support for the life of the product
- Two-year warranty to ensure it remains in top working condition



BMP®21-LAB Printer and Accessories

Part No.	Description
139538	BMP21-LAB Printer
139540	Rechargeable li-ion battery
139542	Hard carry case
875407	Li-ion battery & ac adapter pack
110891	Magnet accessory
110892	Multifunction tool (magnet-light-stand accessory)
110937	AC adapter

Printer ships with: Printer + Lanyard + M21-750-7425 label cartridge

Label Materials that Last

Ultra-thin Polyester (B-461)



- Ultra-thin label easily slides in and out of racks and centrifuges
- Ideal for tubes, vials and conicals
- Lasts in liquid nitrogen, autoclaves and hot water baths

Part No.	Reference	Colour	Width (mm)	Length (m)
110931	M21-375-461	Black on White	9.53	6.40
110932	M21-500-461	Black on White	12.70	6.40
110933	M21-750-461	Black on White	19.05	6.40

Polypropylene Label (B-7425)



- Versatile, durable and economical
- Ideal for vials, tubes and other specimen containers
- Lasts in liquid nitrogen

Part No.	Reference	Colour	Width (mm)	Length (m)
121014	M21-375-7425	Black on White	9.53	6.40
121015	M21-500-7425	Black on White	12.70	6.40
121016	M21-750-7425	Black on White	19.05	6.40

Chemical-resistant Polyester Label (B-488)



- Ideal for slides
- Lasts in freezers and autoclaves, and withstands DMSO, ethanol and short-term xylene exposure

Part No.	Reference	Colour	Width (mm)	Length (m)
110934	M21-375-488	Black on White	9.53	6.40
110935	M21-500-488	Black on White	12.70	6.40
110936	M21-750-488	Black on White	19.05	6.40

Nylon Cloth Label (B-499)



- Easy-to-handle label for non-rigid and curved surfaces
- Ideal for bags, bottles and tubes
- Lasts in liquid nitrogen, freezers and autoclaves

Part No.	Reference	Colour	Width (mm)	Length (m)
110893	M21-375-499	Black on White	9.53	4.88
110894	M21-500-499	Black on White	12.70	4.88
110895	M21-750-499	Black on White	19.05	4.88

2 YEAR
WARRANTY

New Zealand: 0800 1 289 266 www.microanalytix.co.nz

Australia: 1800 643 339 <http://www.microanalytix.com.au>

BR14B21L © 2014 Brady Worldwide Inc. ALL RIGHTS RESERVED

Clear Polyester Material (B-430)



- For flat or large curved surfaces
- Ideal for general purpose bottles, flasks and dishes
- Lasts in freezers and autoclaves

Part No.	Reference	Colour	Width (mm)	Length (m)
139755	M21-250-430*	Black on Clear	6.35	6.40
139746	M21-250-430-WT-CL*	White on Clear	6.35	6.40
110899	M21-375-430	Black on Clear	9.53	6.40
139748	M21-375-430-WT-CL	White on Clear	9.53	6.40
110900	M21-500-430	Black on Clear	12.70	6.40
139747	M21-500-430-WT-CL	White on Clear	12.70	6.40
110901	M21-750-430	Black on Clear	19.05	6.40
139749	M21-750-430-WT-CL	White on Clear	19.05	6.40

* Cartridge is only compatible with the BMP®21-LAB printer.

Durable Vinyl Material (B-595)



- Ideal for equipment, shelves, cabinets and general facility ID
- Available in a variety of colours
- For uneven and rough surfaces

Part No.	Reference	Colour	Width (mm)	Length (m)
139744	M21-250-595-WT*	Black on White	6.35	6.40
142800	M21-375-595-WT	Black on White	9.53	6.40
142807	M21-500-595-WT	Black on White	12.70	6.40
142797	M21-750-595-WT	Black on White	19.05	6.40
139745	M21-250-595-YL*	Black on Yellow	6.35	6.40
142803	M21-375-595-YL	Black on Yellow	9.53	6.40
142799	M21-500-595-YL	Black on Yellow	12.70	6.40
142811	M21-750-595-YL	Black on Yellow	19.05	6.40
142809	M21-375-595-OR	Black on Orange	9.53	6.40
142796	M21-500-595-OR	Black on Orange	12.70	6.40
142804	M21-750-595-OR	Black on Orange	19.05	6.40
139741	M21-375-595-BK	White on Black	9.53	6.40
139742	M21-500-595-BK	White on Black	12.70	6.40
139743	M21-750-595-BK	White on Black	19.05	6.40
142810	M21-375-595-RD	White on Red	9.53	6.40
142795	M21-500-595-RD	White on Red	12.70	6.40
142801	M21-750-595-RD	White on Red	19.05	6.40
142806	M21-375-595-BL	White on Blue	9.53	6.40
142798	M21-500-595-BL	White on Blue	12.70	6.40
142794	M21-750-595-BL	White on Blue	19.05	6.40
139735	M21-375-595-BR	White on Brown	9.53	6.40
139736	M21-500-595-BR	White on Brown	12.70	6.40
139737	M21-750-595-BR	White on Brown	19.05	6.40
142805	M21-500-595-GN	White on Green	9.53	6.40
		White on Green	12.70	6.40
		White on Green	19.05	6.40
139732	M21-375-595-PL	White on Purple	9.53	6.40
139733	M21-500-595-PL	White on Purple	12.70	6.40
139734	M21-750-595-PL	White on Purple	19.05	6.40
139738	M21-375-595-GY	Black on Grey	9.53	6.40
139739	M21-500-595-GY	Black on Grey	12.70	6.40
		Black on Grey	19.05	6.40

* Cartridge is only compatible with the BMP®21-LAB printer.

BRADY